



All superscript numbers (Example: Description ¹) correspond to important notes at the bottom of this page.

Initial Visit/New Patient Consult Code(s): 99241-99245 or 99202-99205 and 93000	Established Patient/Follow Up Code(s): 99211-99215 ¹ and 93000
---	---

Procedure:	CPT Codes:	Description:
EKG ¹	93000	Electrocardiogram

ICD Programming⁶:	
Follow up visit	99211 – 99215
Programming ICD (Single Chamber)	93289, 93282
Programming ICD (Dual Chamber)	93289, 93283
Programming ICD (Multi Chamber/BI-V)	93289, 93284
DFT Testing	93642
Defibrillator	93289
Pacemaker Programming⁶:	
Follow up visit	99211 – 99215
Pacemaker (Single Device)	93288, 93279
Pacemaker (Dual Device)	93288, 93280
Pacemaker BiV (Multi Device)	93288, 93281
Remote Pacemaker Check	93293,93294,93296
Loop Recorder	33285

Important Notes:

1. EKG (93000) to be requested for all visits.